PRISM Stairlifts



PRISM MEDICAL UK

Superior solutions. Exceptional service.

COMFORT. SA

C. Statement

1

FETY. INDEPENDENCE.

At market leading prices and with fast, hassle free installation

Any Prism Stairlift is designed to make your life easier and ensure you continue to enjoy a relaxed, safe and leisurely lifestyle in your own home.

A Prism Stairlift is an excellent product designed to simply provide an easy and convenient way of accessing upper floors quickly and safely, whenever required.

We offer two Prism Stairlift options based on the layout of your staircase – the Prism Horizon Straight Stairlift and the Prism Curve Stairlift. A brand new Prism Stairlift is available to purchase at market leading prices fully fitted with a two year warranty.

Why take needless risks climbing steep stairs when you can use the Prism Stairlift?



PRISM Horizon Straight Stairlift

Safe, reliable and easy to use – at a market leading price and fast installation



The Prism Horizon Straight Stairlift has been developed specifically for straight staircases without curves or bends and will fit unobtrusively and gracefully into your home.

When not in use Prism Straight Stairlifts can be folded away and with a slim profile will not obstruct the stairway for use by other people.

Prism Straight Stairlifts have been ergonomically designed to ensure they are simple and easy to operate.

There are two models to choose from;

- Horizon Standard (maximum user weight capacity of 120kg/19st)
- Horizon Plus (maximum user weight capacity of 140kg/22st with additional features)

Contact Us: To request a quote or assessment, please

Model	Chair	Measurements	Swivel	Footrest	Comments
Horizon Standard	Standard Fixed arm, seat, pad and back Swivel joystick Diagnostic display in carriage Fixed seat height 460mm	MCW: 640mm Space required:	Manual No powered option available Swivel width: 680mm (with arms forward in standard position. Arms can be set back to reduce swivel width.)	Manual No powered option available	Optional powered hinge available Optional outdoor kit available Rail angles 25° - 55°
Horizon Plus	ERGO Adjustable arms, seat, pad and back. Can use Standard or Plus pad. ERGO joystick Diagnostic display in arm Adjustable seat height 465/485/505mm	Weight limit: 140kg (22st) MCW: 640mm Space required: back to knee + 60mm (using ERGO Standard pad. ERGO Plus can save up to 25mm)	Manual Powered option available Swivel width: 680mm (with arms forward in standard position. Arms can be set back to reduce swivel width.)	Lever-linked No powered option available	Optional powered hinge available Optional outdoor kit available Rail angles 25° - 62°

ERGO Seat & Upholstery Options Horizon Plus Only

The height of the user can have a vital impact in the fitting of a stairlift and it's ability to be used safely.

In most cases, the ERGO seat is suitable, however, in cases where space is at a premium, the ERGO Plus seat can help to bring a user's knees further towards their body, reducing the back to knee measurement and potentially making the difference between the user and the stairlift fitting or not.



PRISM Horizon Standard

The Prism Horizon Standard Stairlift delivers all the essential comfort and safety features you need in a neat, slimline package. With a compact profile which sits close to the wall, the Prism Horizon Standard seat comes in neutral colours to blend seamlessly into your home.



Joystick Total control at your fingertips.



Security Lock Prevent unauthorised use of the lift.



Slimline Standard Seat Simple, compact, confortable.



Lever Linked Swivel Turns the user at the top of the stairs.



Rail Angle 25° - 55°

Weight Limit 120kg (19st)



Digital Carriage Display Easy to read, real-time diagnostics.



Lap Safety Belt Additional safety when lift is in operation.



Manual Footrest with Safety Edge Senses obstructions and brings the chair

to a stop.

Contact Us: To request a quote or assessment, please

PRISM Horizon Plus

The Prism Horizon Plus Stairlift is one of the most advanced stairlifts available. It comes with a choice of two seats, the ERGO and the ERGO Plus. The Prism Horizon Plus Stairlift comes with fully adjustable arms, seat pad and back to suit any user and comes with a wealth of built-in features.



Ergonomic Joystick Total control at your fingertips.



Security Lock Prevent unauthorised use of the lift.



ERGO Seat Pad Short or tall, ERGO fits all.



Lever Linked Swivel Turns the user at the top of the stairs. Powered option also available.



Rail Angle 25° - 62°



Weight Limit 140kg (22st)



Digital Arm Display Easy to read, real-time diagnostics.



Retractable Safety Belt Additional safety when lift is in operation.



Lever Linked Footrest with Safety Edge



Adjustable Seat Back Horizontal adjustment of the seat back. Combined with the multi-position seat pad, the Horizon Plus offers multiple configuration options.

PRISM Curve Stairlift

Safe, reliable and easy to use – at a market leading price and fully installed within days



The Prism Curve Stairlift has been developed specifically for a range of curved staircases and will fit unobtrusively and gracefully into your home.

With the revolutionary ergonomically designed chair, a Prism Curve Stairlift can be fully adjusted to ensure the perfect fit for you at installation.

A Prism Curve Stairlift can be easily folded and allows you to swivel the chair on the landing, allowing for safe transfer on and off the stairlift.

There are three models to choose from;

- Curve ERGO (maximum user weight capacity of 120kg/19st)
- Curve ERGO HD (maximum user weight capacity of 160kg/25st) Ideal for plus-size clients
- Curve ERGO SPACE (maximum user weight capacity of 120kg/19st) Ideal for use on narrow staircases

Contact Us: To request a quote or assessment, please

Model	Chair	Measurements	Swivel	Footrest	Comments
Curve ERGO	ERGO Adjustable arms, seat, pad and back. Can use Standard or Plus pad ERGO joystick Diagnostic display in arm Seat height 484mm	Weight limit: 120kg (19st) MCW: 660mm Space required: back to knee + 110mm ext (using ERGO Standard pad. ERGO Plus can save up to 25mm)	Manual Powered option available Swivel width: 680mm (with arms forward in standard position. Arms can be set back to reduce swivel width.)	Manual Lever-linked option available Powered option available	Power swivel and lever-linked footrest cannot be selected together
Curve ERGO HD	ERGO HD Adjustable arms, seat, pad and back. Can use Standard or Plus pad ERGO joystick Diagnostic display in arm Seat height 484mm	Weight limit: 160kg (25st) MCW: 660mm Space required: back to knee + 110mm ext (using ERGO Standard pad. ERGO Plus can save up to 25mm)	Manual No powered option available Swivel width: 680mm (with arms forward in standard position. Arms can be set back to reduce swivel width.)	Manual Lever-linked option available Powered option available	Specially suited for plus-size users Exactly same dimensions as standard ERGO seat - no difference in height or size
Curve ERGO SPACE	ERGO SPACE Adjustable arms, seat, pad and back. Can use Standard or Plus pad ERGO joystick Diagnostic display in arm Seat height 484mm	Weight limit: 120kg (19st) MCW: 660mm Space required: back to knee + 50mm ext & int (using ERGO Standard pad. ERGO Plus can save up to 25mm)	Manual Powered option available Swivel width: 680mm (with arms forward in standard position. Arms can be set back to reduce swivel width.)	Manual Lever-linked option available Powered option available	Power swivel and lever-linked footrest cannot be selected together Saves up to 85mm space

ERGO Seat & Upholstery Options

Available on Curve ERGO, Curve ERGO HD and Curve ERGO SPACE

The height of the user can have a vital impact in the fitting of a stairlift and it's ability to be used safely. In most cases, the ERGO seat is suitable.

However, in cases where space is at a premium, the ERGO Plus seat can help to bring a user's knees further towards their body, reducing the back to knee measurement and potentially making the difference between the user and the stairlift fitting or not. ERGO ERGO Plus

PRISM Curve ERGO

The Prism Curve ERGO Stairlift is one of the most feature rich curve stairlifts available. It is fully adjustable to the size and shape of the user and comes as standard with an ERGO seat. A multitude of upholstery options ensure that the stairlift fits unobtrusively into your home.



Ergonomic Joystick Total control at your fingertips.



Security Lock Prevent unauthorised use of the lift.



ERGO Seat Pad Short or tall, ERGO fits all.



Adjustable Seat Back Horizontal adjustment of the seat back. Combined with the multi-position seat pad, the ERGO offers multiple configuration options.



Weight Limit 120kg (19st)



Digital Arm Display Easy to read, real-time diagnostics.



Retractable Safety Belt Additional safety when lift is in operation.



Lever Linked Footrest with Safety Edge

Open and close with ease. Powered footrest also available.



Safety Edge Senses obstructions and brings the chair to a stop.

Contact Us: To request a quote or assessment, please

PRISM Curve ERGO HD

The Prism Curve ERGO HD offers all the same high quality features of the standard Prism Curve ERGO but with an increased weight capacity of up to 160kg (25st).

Its increased weight capacity makes this stairlift the ideal solution for plus size clients.

Despite the increased weight capacity, size, shape and seat height are identical to the standard model, ensuring that the Curve ERGO HD fits gracefully into your home.







PRISM Curve ERGO Space

The Prism Curve ERGO Space has been specifically designed for those who require a stairlift to be fitted on a narrow staircase. Its unique compact design provides an extra 60mm of space compared to the Prism Curve ERGO.

For those with the narrowest of staircases, use of the ERGO Plus seat pad is advised as this increases the space saving by an extra 25mm to a total of 85mm.

A comfortable ride is a certainty with the Curve ERGO Space as the back rest curves around the user and is 80mm shorter than that of the standard Curve ERGO.







Technical Specifications

Features	Horizon Standard	Horizon Plus	Curve ERGO	Curve ERGO HD	Curve ERGO SPACE
Type of stairlift	Straight	Straight	Curved	Curved	Curved
Choice of upholstery	×	~	 ✓ 	 ✓ 	 ✓
ERGO Standard seat	×	~	v	×	 ✓
ERGO Plus seat	×	~	 ✓ 	×	 ✓
ERGO Space seat	×	×	 ✓ 	×	 ✓
ERGO HD seat	×	×	 ✓ 	 ✓ 	×
Horizon Standard Seat	~	×	×	×	×
Classic seat	×	×	×	×	×
Adjustable chair	×	~	v	×	~
Diagnostic display in arm	×	~	 ✓ 	×	~
Diagnostic display in carriage	~	×	×	×	×
Manual swivel	~	~	 ✓ 	 ✓ 	 ✓
Powered seat swivel	×	~	 ✓ 	×	~
Manual footplate	~	×	 ✓ 	 ✓ 	 ✓
Lever linked footplate	×	~	 ✓ 	 ✓ 	~
Powered footplate	×	×	 ✓ 	 ✓ 	~
Retractable seat belt	×	~	 ✓ 	 ✓ 	~
Lap belt	~	×	×	×	×
Foot covers	×	×	v	 ✓ 	~
Bariatric option	×	×	×	 ✓ 	×
Powered rail hinge	~	~	×	×	×
Outdoor model	~	~	×	×	×

✓ Compliant with the relevant accreditations

2 YEAR WARRANTY Subject to Terms & Conditions

Prism Medical Head Office

Tel: +44 (0) 844 980 2260Fax: +44 (0) 844 980 2261Email: info@prismmedical.co.ukUnit 4, Jubilee Business Park, Jubilee Way, Grange Moor, Wakefield, West Yorkshire, WF4 4TD

Prism Medical Manufacturing Centre

Tel: +44 (0) 844 980 2296Fax: +44 (0) 844 980 2297Email: info@prismmedical.co.ukUnit 1, Tir Llwyd Industrial Estate, St Asaph Avenue, Kinmel Bay, Rhyl, Conwy LL18 5JZ



Unit 4, Jubilee Business Park, Jubilee Way, Grange Moor, Wakefield, West Yorkshire, WF4 4TD Tel: +44 (0) 845 409 2000 E-mail: info@prismmedical.co.uk www.prismmedical.co.uk