

Freeway M205

The Freeway M205 is a powered mobile hoist designed to deliver safe, dignified, and efficient patient transfers in everyday care environments. With a robust lifting capacity of up to 205kg and intuitive controls, it reduces physical strain on caregivers while ensuring comfort and security for the individual being moved. Whether transferring from bed to chair, bath, or floor, the M205 combines reliability, ease of use, and thoughtful design to support high-quality care.



UK Part Number	1301MH-FW-M205 (standard carry bar) 1301MH-FW-M205Q (QRS carry bar)
User weight limit (SWL)	205kg / 32st / 451lbs
Total weight (unboxed)	<40kg
Castor size	75mm (front), 100mm (rear)
Ground clearance	15mm
Under-bed clearance	115mm
Turning circle	1010mm
Length	1235mm
External width	620-1150mm
Internal width	500-1100mm
Height	1340mm-2025mm
Floor to carry bar height	380-1750mm

Reliable emergency functionality

Manual emergency lowering via twist collar ensures users are never stuck in the hoist, even in power failure scenarios.



Smooth mobility across surfaces

Twin castors allow the hoist to move easily over carpeted and hard flooring, enhancing manoeuvrability in varied environments.



Effortless operation

The easy-to-use handset controls raising/lowering and leg opening, with plug-and-play functionality for quick replacement if needed.



Clear and responsive control system

The control box includes an emergency stop, manual raise/lower, and a display screen showing battery status, usage data, and service reminders.



Convenient storage and transport

The hoist can be folded or broken down into two parts without tools, simplifying storage and portability.

Stable and versatile positioning

Powered wide leg opening improves stability and makes transfers to and from large armchairs easier.

Quick accessory changes

The optional QRS Carry Bar allows for fast switching or integration with QRS scales for easy weighing.