

Freeway SA180C

The Freeway SA180C is a compact sit-to-stand aid designed for individuals with limited mobility who can bear partial weight. It offers a gentle, secure lift, and when paired with the Freeway Thoracic Sling, it offers the ultimate sit-to-stand combination. The angled lifting arm reduces underarm pressure, promoting a natural forward-and-upward movement. With support under the feet, legs, and lower back, the SA180C ensures a safe and comfortable transfer experience.

| UK Part Number | 1301SA-SA180C |
|--------------------------------------|-----------------------|
| User weight limit (SWL) | 180kg / 28st / 397lbs |
| Total weight (unboxed) | <52kg |
| Castor size | 75mm |
| Ground clearance | 96mm |
| Under-bed clearance | 132mm |
| Turning circle | 915mm |
| Length | 875mm |
| External width | 570-980mm |
| Internal width | 415-828mm |
| Height | 925-1570mm |
| Floor to pivot support handle height | 869-1490mm |



Moves smoothly over carpeted floors
twin castors with rear brakes enhance manoeuvrability and control



Supports safe sit-to-stand and toileting transfers
compatible with multiple sling types, including thoracic support slings to reduce drag lifting



Accommodates taller users
trombone-style lifting arm provides extended reach



Provides clear usage insights
display screen shows battery status, service reminders, and usage metrics



Ensures emergency readiness
includes manual lowering and emergency stop functions



Simplifies operation
intuitive handset controls boom movement and leg adjustment



Reduces damage risk during transport
control box positioned within the frame for added protection

