

# Freeway Shower Chair configurator

Build and cost up your shower chair. These are our most popular options. If you can't find what you want, please get in touch to discuss your needs.

Step 1

Choose your model

☐ T40 ☐ T40 AUTO ☐ T50 ☐ T60 ☐ T70 ☐ T70 AUTO ☐ Height Adjustable

Step 2

Choose your seat

Select your seat width

☐ A 540mm (21")  
☐ B 490mm (19")  
☐ C 440mm (17")

Select your seat type

☐ Horseshoe

☐ Fully padded

☐ Full ring

Select aperture size

Large

400mm

60mm

210mm

75mm

440mm

☐

Medium

400mm

65mm

160mm

75mm

445mm

☐

Small

400mm

80mm

118mm

65mm

445mm

☐

400mm

60mm

210mm

275mm

65mm

445mm

☐

400mm

65mm

160mm

245mm

90mm

445mm

☐

400mm

80mm

118mm

210mm

110mm

445mm

☐

Select seat extension as required

☐ Extra 50mm ☐ Extra 75mm

Select covering

☐ Dartex ☐ Vinyl

Select seat comfort

☐ Standard ☐ Extra comfort

Please note: All standard seats are fitted with pan sliders.

All AUTO seats are fitted with skirts.

Step 3

Choose your backrest

Select your backrest

Choose between horizontal or vertical push handles, useful when space is restricted or against a WC cistern.

☐ Horizontal push handles

☐ Vertical push handles

Please note: All AUTO chairs come with pram handles as standard.

Step 4

Choose your armrests

☐ Standard

Gives 370mm of clear width on a standard 'C' (440mm) frame.

☐ Half wide

Provide 50mm extra width.

☐ Wide

Provide 100mm extra width.

\* Measurements are taken between armrests across the narrowest point.

Step 5

Choose your footrests

☐ Standard footrests

Adjust from the floor to 195mm.

☐ Raised footrests

60/90/120mm from standard.

☐ Heavy duty footrests

Step 6

Choose the transportability

Many Freeway shower chairs are available as a transportable option. This means the frame can be easily dismantled without tools and is perfect for easy storage or transport.

For example:

T40 Assembled

T40 Disassembled (frame only)

The T50 does not have a transportable option. Please speak to your local representative for a complete list of shower chairs with this option.

Select transportable option?

☐ Yes ☐ No

Step 7

Choose your accessories

A range of accessories is available. Please add the accessories you want here.

Your Freeway Shower Chair is complete!

Please pass your completed form to your local representative or call customer services.

Shower chair configurator | Prism Medical

Prism Medical | Shower chair configurator

# Freeway T80 Shower Chair configurator

Build and cost up your shower chair. These are our most popular options. If you can't find what you want, please get in touch to discuss your needs.

Step 1

Choose your model

☐ T80    ☐ T80 AUTO

Step 2

Choose your seat

Select your seat width

☐ A 540mm (21")  
☐ B 490mm (19")  
☐ C 440mm (17")

Select your seat type

☐ Horseshoe    ☐ Fully padded    ☐ Full ring

Select aperture size

Large

☐

Medium

☐

Small

☐

☐ Large  
  
☐

☐ Medium  
  
☐

☐ Small  
  
☐

Select seat extension as required

☐ Extra 50mm    ☐ Extra 75mm

Select covering

☐ Dartex    ☐ Vinyl

Select seat comfort

☐ Standard    ☐ Extra comfort

Please note: All standard seats are fitted with pan sliders.

All AUTO seats are fitted with skirts.

Step 3

Choose your backrest

The T80 comes with extended backrest and pram bar as standard.

Step 4

Choose your armrests

☐ Standard  
  
Gives 370mm of clear width on a standard 'C' (440mm) frame.

☐ Halfwide  
  
Provide 50mm extra width.

☐ Wide  
  
Provide 100mm extra width.

\* Measurements are taken between armrests across the narrowest point.

Step 5

Footrests

The T80 comes with reclining footrests and heel supports as standard.

Step 6

Choose your accessories

A range of accessories is available. Please add the accessories you want here.

Your Freeway Shower Chair is complete!

Please pass your completed form to your local representative or call customer services.

Get in touch

Telephone: 01924 840100

Email: info@prismmedical.co.uk

Website: www.prismmedical.co.uk

Shower chair configurator | Prism Medical

Prism Medical | Shower chair configurator