

Mackworth Raiser

The Mackworth Raiser is a mobility aid designed to support individuals with limited mobility who retain partial weight-bearing ability. It facilitates safe, active sit-to-stand transfers while encouraging user participation, promoting independence, and helping maintain muscle strength.

Ideal for short-distance movements, such as transfers between beds, chairs, and toilets, the Raiser features ergonomic features including padded, height-adjustable knee supports, multiple grip positions, and an independent braking system. These elements enhance comfort, stability, and control during use.

UK Part Number	1301TA-MWRAIS-1
User weight limit (SWL)	160kg / 25st / 352lbs
Total weight (unboxed)	<20kg
Castor size	125mm braked (mid)
Wheel size	35mm (front and rear)
Ground clearance	21mm
Turning circle	575mm
Length	645mm
Width	592mm
Height	1200m



Supports safe short-distance transfers

designed for movement between beds, chairs, and toilets



Offers flexible grip options

multiple hand holds for both client and carer, including angled grips



Maximises user comfort and stability

padded knee support with adjustable height



Provides secure transfer support

compatible with ladder belts or safety belts



Enhances safety during toilet transfers

non-slip base plate accommodates 17" chair and toilet access

Improves control and manoeuvrability

independent braked castors enable smooth pivoting

Simplifies repositioning for carers

rear kick bar allows easier movement