






Comparing hoists

		Fixed									Removable			
MODEL		Freeway TX Advanced 			Freeway Transactive Xtra 				Mackworth Essentials CT 		Prism CP 		Prism CPP 	
SWL (kg)		130	200	272	130	160	200	270	130	200	130	200	130	200
Track type*	1 - 5	•	•	•	•	•	•	•	•	•	•	•	•	•
	6 - 14				•	•	•	•						
	Prism gantry (FSG/PFS)										•	•	•	•
Charging type	Handset								•	•	•	•	•	•
	In track	•	•	•	•	•	•	•	•	•				
Lifts per charge		130	65	35	120	92	60	33	100	60	100	60	100	60
Carry bar	QRS 3 point (black)	•	•	•					•	•	•	•	•	•
	QRS 1 point (white)								•	•	•	•	•	•
	Standard (black)				•	•	•	•						
Traverse	Manual	•	•	•	•	•	•	•	•	•	•	•	•	•
	Powered	•	•	•	•	•	•	•						
Certifications	TUV	•	•	•					•	•	•	•	•	•
	MDR	•	•	•	•	•	•	•	•	•	•	•	•	•
Options**	QRT	•	•	•	•	•	•	•						
	Powered turntable	•	•	•	•	•	•	•						
	Powered H-system	•	•	•	•	•	•	•						
	Return to charge	•	•	•	•	•	•	•						
	Constant charge	•	•	•										

*For more details on the different track types, please get in touch.

**A single option is available per unit.