

TRANSACTIVEXTRA



130kg, 160kg, 200kg, 270kg versions available



FEATURES

The Transactive Xtra is a safe, reliable and effective ceiling track hoist for carers and users alike. Together with the extensive range of track configurations, the Transactive Xtra can provide the perfect solution to your hoisting requirements.

Available in multiple configurations

Smooth running lift tape Unobtrusive design discreetly blends with surroundings Manual or powered traverse options LCD screen for diagnostics and battery status Manual and electrical emergency lowering Padded carry bar with sling retaining clips for added safety available as standard Ergonomically designed pneumatic hand control, with coloured buttons to indicate direction of movement Manual Lowering & Raising **LCD** Features Lift Counter Battery Level Indicator **Diagnostics** Feature Status Indicator

Emergency Pull Cord





INSTALLATION SOLUTIONS





TRACKS AND ACCESSORIES



The Transactive Xtra is manufactured in the UK and has been tested to ISO 10535

TRANSACTIVE XTRA SPECIFICATIONS

Max. Safe Working Load: 130kg, 160kg, 200kg, 270kg

Fire Retardant ABS Housing

Charging Time: 8 Hours Max.

Unit Weight: Manual - 9.7kg Driven - 10.4kg

Lifting Range: 2000mm

Lifting Speed: Loaded - 30mm/s Unloaded - 56mm/s

SLA Batteries: 2 x 12V, 5Ah

Hoist Protection: IP21 Handset Protection: IP67

Warranty: 36 months (Excluding Batteries)

Prism Medical UK Head Office Unit 4, Jubilee Business Park, Grange Moor, Wakefield WF4 4TD

Telephone: 01924 840100 Email info@prismmedical.co.uk

Prism Medical UK Manufacturing Centre Unit 1, Tir Llwyd Industrial Estate, Kinmel Bay, Conwy LL18 5JZ



