

# Transfer Aids Stand Aids Mobile Hoists



prismmedicaluk  
a prismhealthcare company

# Transfer Aids

A range of transfer aids to meet patient needs for use within the home or care environment manufactured from powered coated steel.



mackworth  
essentials



## Mackworth Essentials Patient Turner

- Compact design
- Transport wheels
- Height adjustable lower leg pad
- Height adjustable handle bars

### Specifications

UK part number	1301TA-MWESSTURN
User weight limit (SWL)	160kg
Total weight (unboxed)	14kg
Castor size	75mm
Length	480mm
Width	430mm
Height (max)	1260mm

mackworth



## Mackworth Raiser

- Base plate shaped to enable the raiser to be manoeuvred around furniture and toilets
- Adjustable height lower leg pads – available in single or dual pad option
- Independent locking brakes on the mid wheel to assist with control and manoeuvrability
- Multiple grip options for both the user and carer
- Can be used with strap and belts

### Specifications

UK part number	1301TA-MWRAIS-1 (single pad) 1301TA-MWRAIS-2 (dual pad)
User weight limit (SWL)	160kg
Total weight (unboxed)	21kg
Wheel size	125mm braked (mid)
Castor size	35mm (front and rear)
Ground clearance	20mm
Length	650mm
Width	590mm
Height	1190mm

freeway



## Freeway TA180 Transfer Aid

- Provides easy transportation around the home
- Allows better eye contact between the user and carer
- Manual leg opening
- Swing away seat pads to enable access and security for the patient
- Positions the patient in an upright position
- Removable foot plate
- Can be used with a safety belt


### Specifications

UK part number	1301TA-FWTA180
User weight limit (SWL)	180kg
Total weight (unboxed)	38kg
Castor type	Braked (rear)
Castor size	100mm
Length	930mm
Width	680mm
Height	1015mm
Turning circle	1010mm
Leg Opening	Manual



# Stand Aids

A range of powered stand aids to meet patient needs to assist in raising to a standing position. Manufactured from powder coated steel and Linak components all with emergency mechanisms the products are robust and provide safety when in use.



mackworth  
essentials



## Mackworth Essentials SA180E

- Robust unit that provides safety when in use
- Manual leg opening
- Removable foot plate
- Soft, adjustable lower leg support with integral strap attachment points

### Specifications

UK part number	1301-ESS180
User weight limit (SWL)	180kg
Total weight (unboxed)	47kg
Lifting range (min-max)	790-1670mm
Castor type	Braked (rear)
Castor size	100mm
Length	1220mm
Width	670mm
Height	1150mm
Leg Opening	Manual

mackworth



## Mackworth SA180

- Robust unit that provides safety when in use
- Powered leg opening and closing
- Removable foot plate
- Soft, adjustable lower leg support with integral strap attachment points
- Battery status alarm

### Specifications

UK part number	1301SA-SA180
User weight limit (SWL)	180kg
Total weight (unboxed)	47kg
Lifting range (min-max)	790-1670mm
Castor type	Braked (rear)
Castor size	100mm
Length	1220mm
Width	670mm
Height	1150mm
Leg Opening	Powered

freeway



## Freeway SA180C

- Compact design
- Powered leg opening and closing
- Removable foot plate
- Soft, adjustable lower leg support with integral strap attachment points
- LCD display with battery status and Service indicator
- Adjustable carry bar with 5 set positions

### Specifications

UK part number	1301SA-SA180C
User weight limit (SWL)	180kg
Total weight (unboxed)	47kg
Lifting range (min-max)	820-1490mm
Castor type	Braked (rear)
Castor size	75mm
Length	880mm
Width	540mm
Height	920mm
Leg Opening	Powered

# Mobile Hoists

A range of mobile hoists to assist in lifting a patient in a home or care environment. Manufactured from Linak components all with emergency mechanisms the products are robust and provide safety when in use.



mackworth  
essentials



## Mackworth Essentials 180 Mobile Hoist

- Compact design
- Powered leg opening
- Powered coated steel frame
- Tool free assembly

### Specifications

UK part number	1301MH-ESS180
User weight limit (SWL)	180kg
Total weight (unboxed)	40kg
Lifting range (min-max)	380-1755mm
Castor type	Braked (rear)
Castor size	75mm front 100mm rear
Length	1210mm
Width	625mm
Height	1330mm
Leg Opening	Powered

mackworth



## Mackworth M180 Fixed/Folding Mobile Hoists

- Available in fixed or folding mast options
- Manual leg opening
- Powered coated steel frame
- Tool free assembly
- Battery status alarm

### Specifications

UK part number	1301MH-M180FIXED 1301MH-M180FOLD
User weight limit (SWL)	180kg
Total weight (unboxed)	39kg
Lifting range (min-max)	220-1660mm
Castor type	Braked (rear)
Castor size	100mm
Length	1310mm
Width	680mm
Height	1360mm
Leg Opening	Manual

freeway



## Freeway S180E Mobile Hoist

- Compact design
- Powered leg opening
- Powered coated steel frame
- Tool free assembly
- LCD display with battery status and Service indicator

### Specifications

UK part number	1301MH-FW-S180E
User weight limit (SWL)	180kg
Total weight (unboxed)	39kg
Lifting range (min-max)	380-1755mm
Castor type	Braked (rear)
Castor size	75mm front 100mm rear
Length	1210mm
Width	625mm
Height	1330mm
Leg Opening	Powered

freeway



## Freeway S205 Mobile Hoist

- Compact design with an increased weight limit
- Powered leg opening
- Powered coated steel frame
- Tool free assembly
- LCD display with battery status and Service indicator

### Specifications

UK part number	1301MH-FW-M205
User weight limit (SWL)	205kg
Total weight (unboxed)	40kg
Lifting range (min-max)	380-1755mm
Castor type	Braked (rear)
Castor size	75mm front 100mm rear
Length	1210mm
Width	625mm
Height	1330mm
Leg Opening	Powered

# Bariatric Hoist

prism

## Prism A320B Mobile Hoist

- 320kg weight limit
- 4 point carry bar
- Coated aluminium frame
- Support bars for gait training
- Powered leg opening
- The mast can be assembled in three different height mounting positions to enhance the maximum or minimum height required
- LCD display with battery status and Service indicator



### Specifications

UK part number	1301MH-PUA-320B
User weight limit (SWL)	318kg
Total weight (unboxed)	64kg
Lifting range (min-max)	492-1785mm
Castor type	Braked (rear)
Castor size	100mm front 125mm rear
Length	1450mm
Width	775mm
Height	1460mm
Leg Opening	Powered

[prismmedical.co.uk](http://prismmedical.co.uk)

Unit 4, Jubilee Business Park,  
Jubilee Way, Grange Moor, Wakefield WF4 4TD

prismmedicaluk  
a prismhealthcare company