

Through Floor Lifts



Prism Medical UK are proud to work in partnership with specialist Through Floor Lift manufacturers, Terry Mobility Solutions and Wessex Medical Equipment.

Through Floor Lifts enable wheelchair users to travel between floors in their wheelchair, avoiding the need to transfer from the wheelchair and onto, for example, a stairlift.



An Introduction to Through Floor Lifts

Regain the freedom of your home by being able to travel between floors smoothly, quietly and safely either seated or in one of our wheelchair carrying models.

The installation of a through floor lift is a far less complicated operation than you might think. Following the initial survey and the subsequent decision to proceed, our experienced team will carry out all necessary preparatory work and install the lift in four to five days providing this is one continuous process.



Lower room when lift upstairs



Upper room when lift downstairs



Entering the lift upstairs



Moving downstairs



Lift approaching lower floor



Exiting the lift downstairs

Key Features

- Powered or manual door options
- Up to 250kg (39 stone) maximum load capacity
- Smooth and quiet operation
- Inherent safety system to guard against free fall of car
- Simple and unintrusive
- Radio wireless controls in car and at both floors which enable travel, stop and door open/close commands
- Adjustable lift car controls can be positioned on the left or right of the car to suit users needs
- No unsightly surface trunking
- A range of optional accessories and customisation

Key Measurements

| Lift Model | Internal Width | Internal Length | External Width (y) | External Length (x) |
|---------------------------------------|----------------|-----------------|--------------------|---------------------|
| Terry Harmony C (Compact/seated) | 680 | 560 | 710 | 915 |
| Terry Harmony S (Std wheelchair) | 1175 | 755 | 920 | 1265 |
| Terry Harmony W (Wider) | 1175 | 890 | 1055 | 1265 |
| Terry Harmony L (Longer) | 1375 | 755 | 920 | 1465 |
| Terry Harmony LW (Longer Wider) | 1375 | 890 | 1055 | 1465 |
| Wessex HomeLift VM30 (Compact/seated) | 590 | 790 | 722 | 1090 |
| Wessex HomeLift VM31 (Std wheelchair) | 740 | 1180 | 872 | 1475 |
| Wessex HomeLift VM36 (Lge wheelchair) | 790 | 1230 | 922 | 1525 |
| Wessex HomeLift VM50 (Compact/seated) | 590 | 790 | 722 | 1090 |
| Wessex HomeLift VM51 (Std wheelchair) | 740 | 1180 | 872 | 1475 |
| Wessex HomeLift VM56 (Lge wheelchair) | 790 | 1230 | 922 | 1525 |

To arrange a demonstration of these Through Floor Lifts please contact us on 0845 409 2000 or visit www.prismmedical.co.uk