

PRISM GK Assisted Bathing Systems



PRISM MEDICAL UK

Moving, Handling and Bathing Solutions™



SUPERIOR A



ASSISTED BATHING SYSTEMS.

Highest quality build at value for money prices

The Prism GK Assisted Bath Range provides a multitude of bathing solutions for a wide variety of environments, including hospitals, care homes and individual residences.

Safety, quality and reliability have been at the forefront of the design, ensuring that your clients benefit from only the finest bathing system.

Included within the range are four models; Prism GK Relax, Prism GK Oasis, Prism GK Magic 2000 and the Prism GK Compact.

Each bath has been developed in partnership with GK, a long established Austrian manufacturer. The Prism GK Assisted Bathing Systems have been specifically designed to provide a responsive

solution to cater for a number of clients, from low to very high levels of dependency.

All baths comply with WRAS (Water Regulatory Advisory Scheme) standards, so you can be assured that they are manufactured to the highest quality in accordance with current industry regulations.



Available as added extras on the Prism GK Assisted Bathing Systems:

Hydro-massage therapy

Hydro-massage therapy is an ultra-safe usage of air spa stimulation, designed to aid rehabilitation processes, promote good blood circulation and relax/invigorate the sensory system.

Colour Therapy

Colour therapy is a dynamic system utilising over 100 coloured LED's, designed to colour the bathwater and provide a relaxing and atmospheric bathing experience.

PRISM GK Oasis

Prism GK Oasis Bath – with integrated hoisting system.

The Prism GK Oasis bath offers a host of innovative features to ensure a comfortable bathing environment for the client and to provide a safe working area for the carer.

This state of the art bath accommodates a full range of dependency levels with its variable height and integrated hoist, ensuring transfers into and out of the bath can be completed with ease.



To request a quote or assessment, call us on 0845 409 2000 or visit our website www.prismmedical.co.uk

The Prism GK Oasis bath is suitable for nursing and residential homes.

SAFE WORKING LOAD OF HOIST
150KG

Features and benefits



Touch free control panel

State of the art touch free control panel minimises cross infection and ensures optimal bathing hygiene.

Integrated hoisting system

Eliminates the need for a separate hoist and creates more usable space around the bath.

Variable height adjustment

Enabling the bath to be lowered and raised at the touch of a button depending on the needs of both the client and care giver.

Optional hoist chair

Works seamlessly alongside the bathing system for use when locating the client from another room.

Hydro-massage therapy

The Prism GK Oasis can be fitted with a spa style Hydro-massage therapy system to provide a first class bathing experience.



Ergonomic controls

Allows easy usage of all bathing features by the care giver.

Safe water temperature control

Water temperature controlled with TMV 3 approved mixer valves to prevent scalding.

Integrated shower attachment

Suitable for hair and body washing, proving the client with the most cleansing bathing experience.

Colour therapy

Optional colour therapy treatment available to provide a unique relaxing and atmospheric bathing experience.

BATH TUB MAXIMUM LIFTING CAPACITY
(Resident + water) 400KG

DEPENDENCY LEVEL:
LOW TO MEDIUM

EXTENSIVE 5 YEAR WARRANTY AVAILABLE

PRISM GK Relax

Prism GK Relax – with vertical raise entrance door for easy access.

The Prism GK Relax recline bath is designed to be used in a multitude of bathing areas. The Prism GK Relax bath is suitable for nursing and residential homes.

Ease of access is the key behind this particular design, as it incorporates an intuitive side entry system to enable seamless transfers for the client.

The Prism GK Relax bath has an auto fill function which allows a set volume of water to fill the bath at a desired temperature, saving time which can be dedicated to the all-important patient care.

To request a quote or assessment, call us on 0845 409 2000 or visit our website www.prismmedical.co.uk



The Prism GK Oasis bath is suitable for nursing and residential homes.

Features and benefits



Vertical raise entrance door

Provides easy access either independently, or utilising mobile or ceiling track hoists.

Reclining system

Works in tandem with the in built moulded seat, making this the ideal bathing solution for clients with poor trunk or postural control.

Audio/sound therapy

The Prism GK Relax can be fitted with an audio/sound therapy system which uses music as a highly beneficial aid in the treatment process.

Hydro-massage therapy

The Prism GK Relax can be fitted with a spa style Hydro-massage therapy system to provide a first class bathing experience.



Touch screen control panel

Centralised control panel allows the bath to be seamlessly operated by either the client or care giver.

Safe water temperature control

Water temperature controlled with TMV 3 approved mixer valves to prevent scalding.

Variable height adjustment

Enabling the bath to be lowered and raised at the touch of a button depending on the needs of both the client and care giver.

Colour therapy

Optional colour therapy treatment available to provide a unique relaxing and atmospheric bathing experience.



BATH TUB MAXIMUM LIFTING CAPACITY
(Resident + water) 400KG

DEPENDENCY LEVEL:
LOW TO VERY HIGH

EXTENSIVE 5 YEAR WARRANTY AVAILABLE

PRISM GK Magic 2000

Prism GK Magic 2000 – with touch screen control panel.

The Prism GK Magic 2000 is a highly flexible height adjustable bath incorporating the latest state of the art technology.

The ergonomically designed bath ensures greater movement for the user, with an option of two bath lengths to provide that bespoke feel.


Ultimate comfort can be achieved using a multitude of in-built settings available via the touch screen control panel.

To request a quote or assessment, call us on 0845 409 2000 or visit our website www.prismmedical.co.uk



The Prism GK Magic 2000 is suitable for acute hospital environments and nursing and residential homes.

Features and benefits



Touch screen control panel

Centralised control panel provides easy access to all bathing functions.

Automatic fill

The Prisma GK Magic 2000 is fitted with 3 automatic fill levels, allowing the client to relax while the bath fills effortlessly.

Audio/sound therapy

The Prisma GK Magic 2000 can be fitted with an audio/sound therapy system which uses music as a highly beneficial aid in the treatment process.

Hydro-massage therapy

The Prisma GK Magic 2000 can be fitted with a spa style hydro-massage therapy system to provide a first class bathing experience.

Integrated shower attachment

Suitable for hair and body washing, proving the client with the most cleansing bathing experience.

Multiple length options

The Prisma GK Magic 2000 is available in two lengths, ensuring the perfect fit for the client.

Variable height adjustment

Enabling the bath to be lowered and raised at the touch of a button depending on the needs of both the client and care giver.

Colour therapy

Optional colour therapy treatment available to provide a unique relaxing and atmospheric bathing experience.

BATH TUB MAXIMUM LIFTING CAPACITY
(Resident + water) 400KG

DEPENDENCY LEVEL:
MEDIUM TO VERY HIGH

EXTENSIVE 5 YEAR WARRANTY AVAILABLE

PRISM GK Compact

Prism GK Compact – with easy to operate control panel.

The Prism GK Compact height adjustable bath provides an ideal solution for a multitude of dependency needs.

The variable height adjustment allows optimal positioning of mobile hoists to aid client transfers in and out of the bath, whilst the slide in foot board ensures that the bath can be tailored to suit the height of any client.



To request a quote or assessment, call us on 0845 409 2000 or visit our website www.prismmedical.co.uk

The Prism GK Compact bath is suitable for nursing and residential homes.

Features and benefits



Variable height adjustment

Allowing the bath to be raised to a comfortable height for the carer and lowered to aid client transfers.

Slide in foot board

To adjust the length of the bath to suit the height of the client.



Bath fill temperature display

To ensure the patient baths in the optimal required temperature.

Integrated arm rests

For increased client comfort during use.

Easy to operate control panel

Allows for centralised control of the bath. Easy to use controls aid the bathing experience.

Adjustable bath feet

Allows the bath to be levelled and correctly positioned on installation.

Access under the bath

Allowing the positioning of a mobile hoist to aid client transfers.

BATH TUB MAXIMUM LIFTING
CAPACITY
(Resident + water) 400KG

DEPENDENCY LEVEL:
LOW TO MEDIUM

**EXTENSIVE 5 YEAR
WARRANTY AVAILABLE**

Prism GK Bath Range Summary

	Compact	Oasis	Magic 2000	Relax
				
Bath Length (L)	207 cm	172 cm	224 cm	183 / 205 cm
Bath Width (W)	82 cm	79 cm	93 cm	78 cm
Variable Height	●	●	●	●
Tilt Function				●
Patient Shower	●	●	●	●
Integral Hoist		●		
Auto-stop (when filling)	Optional	●	●	●
Length Options			●	
Hydro-massage therapy	Optional	Optional	Optional	Optional
Colour Therapy	Optional	Optional	Optional	Optional
Audio/Sound Therapy	Optional	Optional	Optional	Optional
WRAS Approved	●	●	●	●
TMV3 Approved Mixer	●	●	●	●
Bath tub maximum lifting capacity (resident + water)	400 kg	400 kg	400 kg	400 kg
Safe Working Load (SWL) of hoist		150 kg		
Emergency lowering	Optional	●	●	emergency open door / tilt forward

✓ Compliant with the relevant accreditations

5 YEAR WARRANTY AVAILABLE

Subject to Terms & Conditions



Suitable applications/environments

Prism GK Compact – Nursing/residential

Prism GK Oasis – Nursing/residential

Prism GK Magic 2000 – Acute & nursing/residential

Prism GK Relax – Nursing/residential

Your Prism Medical regional office...

Westholme Division

Unit 4, Jubilee Business Park, Jubilee Way, Grange Moor, Wakefield, West Yorkshire
WF4 4TD UK
Tel: +44 (0) 844 980 2260
Fax: +44 (0) 844 980 2261
E-mail: info@westholme.co.uk

Movement2 Division

Unit 8 The Warren, East Goscote Industrial Estate, East Goscote, Leicestershire LE7 3XA UK
Tel: +44 (0) 844 980 2281
Fax: +44 (0) 844 980 2282
E-mail: mail@movement2.co.uk
www.movement2.co.uk

Carefree Division

13 - 15 Pennington Street, Hindley, Wigan, Lancashire WN2 3AH UK
Tel: +44 (0) 844 980 2266
Fax: +44 (0) 844 980 2267
E-mail: info@carefreelifting.co.uk

Prism South Division

Unit 42 South Hampshire Industrial Park, Totton, Southampton, Hampshire SO40 3SA UK
Tel: +44 (0) 844 980 2272
Fax: +44 (0) 844 980 2273
E-mail: prismsouth@prismmedical.co.uk



PRISM MEDICAL UK



MIX Paper from responsible sources
FSC® C003919



Unit 1, Tir Llwyd Industrial Estate, St Asaph Avenue, Kinnel Bay, Rhyl, Conwy, LL18 5JZ UK
Tel: +44 (0) 845 409 2000 E-mail: info@prismmedical.co.uk www.prismmedical.co.uk