

# FREEWAY Raiser



PRISM MEDICAL UK

*Superior solutions. Exceptional service.*

# FREEWAY Raiser

The **Freeway Raiser** is a sit to stand platform, which is designed to provide an active, safe and comfortable transfer of the user. It is ideal for short distance transfers, such as to/from bed to chair.

Its versatility allows use with wheelchairs, toilet and shower chairs.

Soft and wide padding in the leg support improves stability, and the sturdy handle offers many grip options for both the user and the carer.

The rear wheels with independent locking manual brakes, makes it easy for a carer to manoeuvre, and provides maximum stability and safety for the user.

**The Freeway Raiser has a safe working load of 160 kg.**

To meet the needs of the user, two different styles of safety belts are available – **Freeway Raiser Ladder Strap** and **Freeway Safety Belt**. Both can be easily attached to the transfer platform to enable either one or two carers to assist a user to their feet.

The unique handle mechanism on the Freeway Raiser provides further security and safety for belt fastening.



Manufactured in the UK

## Key Benefits and Features

- Helps to assist a user with sit to stand transfers
- Maximum user weight: 160 kg (25 stone)
- Lightweight and easy to manoeuvre
- Padded leg support for added comfort
- Multiple grip options for both user and carer
- Raiser Strap available to enable a carer to assist a user to their feet
- Independent locking rear brakes

## Technical details:

Maximum weight capacity	160 kg / 352.7 lbs
Turning radius	360°
Overall height	1201 mm / 47.3"
Ground clearance	19.5 mm / 0.8"
Front and rear castors	35 mm / 1.4"
Wheels	125 mm / 4.9"
Base assembly	12.5 kg / 27.6 lbs
Assembled unit	19 kg / 41.9 lbs

**For more information please call our team on  
01924 840100 or visit [www.prismmedical.co.uk](http://www.prismmedical.co.uk)**



**PRISM MEDICAL UK**