






Comparing hoists

| | | Fixed | | | | | | | | | Removable | | | |
|------------------|------------------------|---|-----|-----|---|---------------------|---------------------|---------------------|---|-----|---|-----|---|-----|
| Model | | Freeway TX Advanced | | | Freeway Transactive Xtra | | | | Mackworth Essentials CT | | Prism CP | | Prism CPP | |
| | |  | | |  | | | |  | |  | |  | |
| SWL (kg) | | 130 | 200 | 272 | 130 | 160 | 200 | 270 | 130 | 200 | 130 | 200 | 130 | 200 |
| Track type* | 1 - 5 | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | 6 - 14 | | | | • | • | • | • | | | | | | |
| | Prism gantry (FSG/PFS) | | | | | | | | | | • | • | • | • |
| Charging type | Handset | | | | | | | | • | • | • | • | • | • |
| | In track | • | • | • | • | • | • | • | • | • | | | | |
| Lifts per charge | | 130 | 65 | 35 | 120 | 92 | 60 | 33 | 100 | 60 | 100 | 60 | 100 | 60 |
| Carry bar | QRS 3 point (black) | • | • | • | | | | | • | • | • | • | • | • |
| | QRS 1 point (white) | | | | | | | | • | • | • | • | • | • |
| | Standard (black) | | | | • | • | • | • | | | | | | |
| Traverse | Manual | • | • | • | • | • | • | • | • | • | • | • | • | • |
| | Powered | • | • | • | • | • | • | • | | | | | | |
| Certifications | TUV | • | • | • | | | | | • | • | • | • | • | • |
| | MDR | • | • | • | • | • | • | • | • | • | • | • | • | • |
| Options** | QRT | • | • | • | • | • | • | • | | | | | | |
| | Powered turntable | • | • | • | • | • | • | • | | | | | | |
| | Powered H-system | • | • | • | • | • | • | • | | | | | | |
| | Return to charge | • | • | • | • | • | • | • | | | | | | |
| | Constant charge | • | • | • | • (track 1 only) | • (track 1 only) | • (track 1 only) | • (track 1 only) | | | | | | |

* For more details on the different track types, please get in touch.

**A single option is available per unit.