





























Technical Specifications















Below are the technical specifications for the most popular Freeway Shower Chair models. Please consult our website, www.prismmedical.co.uk, for the full range of dimensions.










Standard Chairs


These chairs have been weight tested to 200kg (30 stone). Excludes the T100, which is tested to 350kg (55 stone).


	1 - Overall maximum height		8 - Standard depth of the seat
	2 - Overall maximum width		9 - Top of the armrest to top of the seat
	3 - Length with footrests fitted at the lowest position		10 - Width clearance
	4 - Length with footrests fitted at the highest position		11 - Top of backrest handle to top of the seat
	5 - Width of standard seat		11a - Top of backrest vinyl to top of the seat
	6 - Height from floor to top of the seat		12 - Distance between armrest at the rear
	7 - Height clearance		13 - Distance between armrest at the front


T30														
T30A	1025	660	570	570	540	480	355	400	190	485	490	450	475	520
T30B	975	610	570	570	490	480	355	400	190	435	490	450	425	470
T30C	925	560	570	570	440	480	355	400	190	385	490	450	375	420

T40														
T40A	980	640	1050	950	540	560	445	400	190	485	490	450	475	520
T40B	980	610	1050	950	490	560	445	400	190	435	490	450	425	470
T40C	980	560	1050	950	440	560	445	400	190	385	490	450	375	420
T40A AUTO	980	640	1050	950	540	560	480	400	190	485	490	450	475	520
T40B AUTO	980	610	1050	950	490	560	480	400	190	435	490	450	430	470
T40C AUTO	980	560	1050	950	440	560	480	400	190	385	490	450	410	410
T40A COM	980	640	1050	950	540	560	470	400	190	485	490	450	475	520
T40B COM	980	610	1050	950	490	560	470	400	190	435	490	450	425	470
T40C COM	980	560	1050	950	440	560	470	400	190	385	490	450	375	420


T50														
T50A	980	730	960	930	540	560	445	400	190	485	490	450	475	520
T50B	980	680	960	930	490	560	445	400	190	435	490	450	425	470
T50C	980	630	960	930	440	560	445	400	190	385	490	450	375	420

T60														
T60A	980	690	1060	960	540	560	445	400	190	485	490	450	475	520
T60B	980	640	1060	960	490	560	445	400	190	435	490	450	425	470
T60C	980	590	1060	960	440	560	445	400	190	385	490	450	375	420

T70														
	1	2	3	4	5	6	7	8	9	10	11	11a	12	13
T70A	980	750	1100	1150	540	560	445	400	190	485	490	450	475	520
T70B	980	700	1100	1150	490	560	445	400	190	435	490	450	425	470
T70C	980	650	1100	1150	440	560	445	400	190	385	490	450	375	420
T70A AUTO	980	790	1060	1140	540	560	480	400	190	485	490	450	475	520
T70B AUTO	980	740	1060	1140	490	560	480	400	190	435	490	450	430	470
T70C AUTO	980	690	1060	1140	440	560	480	400	190	385	490	450	410	410
T70A COM	980	775	1020	1110	540	560	470	400	190	485	490	450	475	520
T70B COM	980	725	1020	1110	490	560	470	400	190	435	490	450	425	470
T70C COM	980	675	1020	1110	440	560	470	400	190	385	490	450	375	420

T90 Paediatric														
	1	2	3	4	5	6	7	8	9	10	11	11a	12	13
T90 Assistant-propelled	910	620	690	N/A	360	440	N/A	330	140	N/A	480	440	300	340
T90 Self-propelled	910	620	740	N/A	360	440	N/A	330	140	N/A	480	440	300	340
T90WC Assistant-propelled	980	530	740	N/A	360	560	440	330	160	440	480	440	320	325
T90WC Self-propelled	980	690	870	N/A	360	560	440	330	160	440	480	440	320	325

T100 Plus-Size														
	1	2	3	4	5	6	7	8	9	10	11	11a	12	13
T100	970	860	1220	1170	790	560	440	550	190	740	480	440	720	750
T100FR	970	910	1220	1170	790	560	440	550	190	740	480	440	720	750

Height Adjustable														
	1	2	3	4	5	6	7	8	9	10	11	11a	12	13
Height Adjustable A	980	640	1050	970	540	560	445	400	190	485	490	450	475	520
Height Adjustable B	980	610	1050	970	490	560	445	400	190	435	490	450	425	470
Height Adjustable C	980	560	1050	970	440	560	445	400	190	385	490	450	375	420

This frame adjusts in height. The floor to top of seat measurement can be 510, 535 or 560mm.

T80 Chairs

The T80 has been weight tested to 160kg (25 stone).

- 1** - Height
- 2** - Height with detachable headrest
- 3** - Width
- 4A** - Length in an upright position
- 4B** - Length with max. recline, including footrests at the highest position
- 4C** - Length max. recline and detachable headrest fitted
- 5** - Seat width
- 6** - Height to the top of the seat
- 7** - Height clearance
- 8** - Seat depth (standard seat)
- 9** - Top of the armrest to top of the seat
- 10** - Width clearance
- 11** - Top of the backrest to top of the seat (without detachable headrest fitted – with headrest add 300mm)
- 12** - Distance between armrests at rear
- 13** - Distance between armrests

T80															
	1	2	3	4A	4B	4C	5	6	7	8	9	10	11	12	13
T80A	1140	1460	650	900	1150	1300	540	570	445	400	190	480	535	475	520
T80B	1140	1460	600	900	1150	1300	490	570	445	400	190	430	535	425	470
T80C	1140	1460	550	900	1150	1300	440	570	445	400	190	380	535	375	420
T80A COM	1140	1460	730	900	1150	1300	540	570	470	400	190	480	535	475	520
T80B COM	1140	1460	680	900	1150	1300	490	570	470	400	190	430	535	430	470
T80C COM	1140	1460	630	900	1150	1300	440	570	470	400	190	380	535	410	410
T80A AUTO	1140	1460	730	900	1150	1300	540	570	470	400	190	480	535	475	520
T80B AUTO	1140	1460	680	900	1150	1300	490	570	470	400	190	430	535	425	470
T80C AUTO	1140	1460	630	900	1150	1300	440	570	470	400	190	380	535	375	420
T80P	1024	1320	560	760	1060	1220	360	570	470	320	160	500	510	400	400

Freeway Flow

The cradle has been weight tested to 130kg (20 stone).

- A** - Seat height
- B** - Cradle end height
- C** - Handle height
- D** - Overall height
- E** - Base frame length
- F** - Back length
- G** - Seat length
- H** - Leg length
- I** - End length
- J** - Internal width
- K** - Overall width
- L** - Upright length
- M** - Reclined length

Freeway Flow	A*	B	C*	D*	E	F	G	H	I	J	K	L	M
Small	592	176	940	1106	800	620	300	360	150	440	502	950	1260
Medium	592	116	1050	1214	800	750	350	450	150	440	502	1134	1508
Large	592	116	1050	1214	800	750	350	450	150	520	582	1134	1508

All measurements in mm.
*Measurements taken at lowest height setting.