

Freeway Bed Management System

The Freeway bed management system provides an easy-to-use, low-friction solution for patients who require frequent repositioning.



 **MANUFACTURED
IN THE UK**

 01924 840100


 info@prismmedical.co.uk

 www.prismmedical.co.uk

Freeway Bed Management System

The Freeway bed management system consists of two separate units: a nylon, water-resistant undersheet with wide, non-slip sides and a fully breathable, water-proof, quilted polyester softshell top sheet.

The system, designed to provide comfort for patients, also significantly reduces the incidence of pressure ulcers. It is easy for carers to use, reducing their stress and ensuring the best care for the patient. It forms an ideal bed solution as part of a single-handed care package.



Features and Benefits

- Provides comfort for patients
- Minimal effort from carers
- Ideal for use by a single carer
- Works easily with four-section profiling beds
- Can be used on a pressure care mattress
- Can also be used with a hoist to assist with turning
- Dimensions: 140 x 200cm
- Maximum user weight: 200kg for standard, 300kg for bariatric

Freeway Positioning Wedge

The Freeway positioning wedge, a companion to the Freeway bed management system, offers firm support for patients in bed. This feature reassures carers and healthcare professionals, ensuring that the patient remains in the required position, promoting their comfort and well-being.



Features and Benefits

- Provides a positive, non-slip action when used against cotton bed linen
- Can be used to help reduce the incidence of pressure ulcers
- Made from 100% soft and flame-retardant foam
- The polyurethane cover is machine washable up to 80°C
- The nylon cover can be wiped clean with disinfectant
- Available in standard and bariatric sizes

Please get in touch to learn more about our moving and handling solutions.

	Product Code
Freeway bed management system (complete)	2002ACCBMS-COM
Freeway bed management system (base sheet only)	2002ACCBMS-BSS
Freeway bed management system (quilted top sheet only)	2002ACCBMS-QS
Freeway bed management system - bariatric (complete)	2002ACCBMS-BCOM
Freeway bed management system - bariatric (base sheet only)	2002ACCBMS-BSB
Freeway bed management system, - bariatric (quilted top sheet only)	2002ACCBMS-QB
Freeway positioning wedge (53cm wide)	2002FW-WED53
Freeway bariatric positioning wedge (60cm wide)	2002FW-WED60B

Prism Medical UK Head Office

Unit 4, Jubilee Business Park
Grange Moor, Wakefield
WF4 4TD

Telephone: **01924 840100**

Email: **info@prismmedical.co.uk**

Website: **www.prismmedical.co.uk**

